# SPL REFERENCE MATERIALS

# CERTIFICATE Alloy Steel QCM for Spectrometry SP 1 – SP 8

#### Status

The Quality Control Materials – QCM comply with the latest ISO Guide 35 definition of the Reference Material.

### Purpose

The QCM are primarily intended for quality assurance purposes in the routine spectrometric analysis of special (ie.beyond main-stream) alloy steels, that is to control the state of statistic regulation of continuously operating analysers ("setting-up") to verify the matrix-match in metal sorting etc. They may not substitute the CRM in establishing traceability of values.

# Source and Manufacture

The candidate materials were either cast to ingots with head, bottom and centre removed and the rest machined to final diameter, or meticulously selected from the commercially available rolled steel bars.

#### Supply units

The QCM of 35 to 43 mm in diameter are regularly supplied in heights of 10 and 25mm (as set or individually), with the option of longer samples upon request.

## Homogeneity and Stability

Both random and trend (axial, radial) inhomogeneity was tested in compliance with ISO Guide 35 by various spectrometric techniques of excellent repeatability, and found statistically insignificant. The QCM are stable by the nature of material.

### Characterisation

The QCM were characterised by an interlaboratory experiment with participation of the most competent and experienced industrial and application laboratories. Various spectrometric methods were applied for each value (AES with spark and glow discharge excitation, XRF), supported when necessary by alternative techniques (wet-way, combustion, thermoevolution). The **participating laboratories** were:

Leco Instrumente Plzeň, Plzeň Mittal Steel Ostrava, Ostrava PCS, Žďár nad Sázavou PSP Slévárna, Přerov Škoda Kovárny, Plzeň

Thermo ARL, Ecublens Třinecké Železárny, Třinec Vítkovice Testing Center, Ostrava Západočeská Univerzita v Plzni, Plzeň Žďas, Žďár nad Sázavou

### Assigned values

are robust means of at least five laboratory means accepted by the technical assessment. The values of the succesive batches, distinguished by the last capital of the QCM code are traced to the original batch values by three qualified laboratories. The particular **uncertainties** are not stated, yet their estimates do not exceed double

of those for the corresponding elements/contents in the contemporary matrixmatching CRM. They are reflected by rounding of the assigned values.

The assigned values are **traceable** to the adequate CRM (ČKD, BAS, Brammer Standards and other) only, no direct traceability has beed established.

QCM	C	Mn	Si	Р	S	Cu	Cr	NI	ΔΙ	84
SP-1A	0.047	1 87	0.33	0.024	0.26	0.52	47.7		AI	IVIO
SP-3B	0.27	0.20	0.70	0.024	0.20	0.52	17.7	8.6	0.004	0.42
CD 4D	0.21	0.29	0.72	0.023	0.008	0.62	15.1	5.65	0.08	0.24
SP-4B	0.31	1.36	1.51	0.027	0.008	0.04	20.8	35.8	0.025	0.02
SP-5B	0.20	1.86	3.07	0.108	0.023	0.15	0.38	3 00	0.18	0.02
SP-6A	C.10	0.38	4.65	0.017	0.009	0.12	0.11	0.04	0.32	0.10
SP-7A	0.006	0.08	0.036	0.007	0.010	0.08	0.01	173	0.02	0.01
SP-8B	2.37	0.86	1 40	0.022	0.012	0.075	27.0	47.5	0.003	0.01
OCM	14/	0.00	1.40	0.022	0.012	0.075	37.0	2.72	0.13	0.10
QUIVI	VV	V		Co	As	Sn	В	Nb	Pb	Sb
SP-1A	0.03	0.058	0.02	0.095	0.006	0.01	0.0007	0.012		
SP-3B	0.12	0.10	0.13	0.02		0.01	0.88			
SP-5B	0.62	0.71	0.35	0.135	0.19	0.08	0.14	0.09	0.09	0.07
SP-6A	0.02	0.016	0.008	0.003	0.003	0.01	0.11	0.00	0.00	0.07
SP-7A	7	0.001	0.004	0.003						
SP-8B	0.05	0.13	0.13	0.075	0.05	0.06	0.03	0.04		

# Assigned values in % m/m<sup>\*</sup>

\*sum of non-listed TMI below 0.05%

# **Users instruction**

The working surface of the QCM must be prepared before the analysis in accordance with the particular analyser manual. The storage in dry and non-corrosive environment is recommended. There are no safety hazzards in the storage and proper use of QCM.

### Producer

SPL is the authorised producer of the CRM for the Czech Metrology Institute, producer of its own RM, QCM and laboratory consumables, organiser of the Proficiency Testing Programs on behalf of the Czech Accreditation Institute, expert in sample management and provider of sampling devices, with the QM compliant to ISO 17025 and ISO Guide 34.

Address: SPL, ul. 1. máje 432, CZ - 735 31 Bohumín, Czech Republic Tel./fax +420 596014627, e-mail: <u>info@spl-bohumin.cz</u>

## Responsible person

Star Bohumín, CZ

Vladimír Bogumský, SPL manager, member of the Czech Committee for Chemical Metrology

Issued in Bohumín on March 22, 2013

liby pro laboratoře

CO: 46605134 DIČ: CZ515414033

Valid till March 22, 2023